





规格型号	A/mm	B/mm	C/mm	D/mm	E/mm	F/mm	G/mm
HGZT5015D	5000	1500	800	300	950	1200	200
HGZT6015D	6000	1500	800	300	1150	1200	200

DRAWING PROJECTION 		OMTOOLS	
REV	A0	HGZT系列阻尼平台 (12支腿)	
DRAWN	TODD		
DATE	2022/3/19	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		kg	1:30