






M3 MOUNTING HOLES
THROUGH 12 PLACES


95mm
[3.74in]

SPRING CLIP TO SECURE
PETRI DISHES AND GLASS SLIDES
2 PLACES



45mm
[1.77in]


ROTATION LOCKING SCREW
2 PLACES



50mm[1.97in] TRAVEL MICROMETER
2 PLACES

M6 MOUNTING HOLES
THROUGH 32 PLACES

DRAWING
PROJECTION



OMTOOLS

REV A0

DRAWN TODD

DATE 2024/1/8

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

OMYE200Z

WEIGHT

3.08kg

SCALE

1:4