




OPEN-ARM DESIGN ALLOWS FOR EASY
DROP-IN INSTALLATION & REMOVAL
OF LENS TUBE FROM OPTICAL SETUP



COUNTERBORE FOR M6
FOR CONNECTORS WITH M6 THREADS
AND M4 CAP SCREW



CAPTIVE LOCKING SCREW
FOR SECURING LENS TUBE
2mm HEX



DRAWING PROJECTION		 	OMTOOLS	
REV	A0	OMJC31K		
DRAWN	TODD			
DATE	2023/7/4	WEIGHT	SCALE	
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.02kg	1:1	