



M4 THERAD
12.5mm [0.49in] MA X DEPTH
16 PLACES

CLEARANCE BORE FOR USE WITH
M6 [1/4-20] CAP SCREW
4 PLACES
ALIGNMENT HOLE FOR USE
WITH Ø3.0mm DOWEL PINS
FOR CENTERING OF SLIDES
3 PLACSE [PINS NOT INCLUDED]

M3 THERAD
4 PLACES

M6 THERAD
12.5mm [0.49in] MA X DEPTH
40 PLACES

SM1(Ø1.035-40mm) SERIES
EXTERNAL THREAD
OMSPJ25A RETAINING RING INCLUDED
2.0MM [0.08in] THICK
Ø23mm[0.91in] CLEAR APERTURE

13mm[1/2in] TRAVEL MICROMETER
2 PLACES

DRAWING PROJECTION		OMTOOLS	
REV	A0		
DRAWN	HSY	WEIGHT	SCALE
DATE	2021/8/5	0.70kg	1:2
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES			