


THRU TO CENTER SM2(2.035in-40) SERIES  
INTERNAL THREAD ACCEPTS 50.8mm[2.0in]  
MAXIMUM 7.7mm[0.3in] THICK OPTICS  
2 GYKH50A RETAINING RING INCLUDED  
48.3mm [1.9in] CLEAR APERTURE



A-A

DRAWING PROJECTION

**OMTOOLS**

REV A0

F2H50B

DRAWN TODD

DATE 2020/5/13

WEIGHT

SCALE

FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

0.03kg

1:1